

GNR ANALYTICAL INSTRUMENTS



SUPPLY · SERVICE · SOFTWARE

GNR ANALYTICAL INSTRUMENTS UK - web: www.gnranalytical.co.uk - email: info@gnranalytical.com 0115 937 6686

	Ni 01
	Ni
	Orientation
Al	0.005-7
B	0.005-3
C	0.005-1
Co	0.005-25
Cr	0.005-50.0
Cu	0.005-35
Fe	0.005-50
Mg	0.005-0.1
Mn	0.005-8
Mo	0.005-35
Nb	0.005-6
Ni	Ref.
Si	0.005-7
Sn	0.005-1
Ti	0.005-6
V	0.005-1
W	0.005-12

LIMITS OF DETECTION FOR THE GNR MINILAB 150 in NICKEL BASED MATERIALS