

GNR ANALYTICAL INSTRUMENTS



SUPPLY · SERVICE · SOFTWARE

GNR ANALYTICAL INSTRUMENTS UK - web: www.gnranalytical.co.uk - email: info@gnranalytical.com 0115 937 6686

MG 01 MINILAB 300 Mg Orientation	
Ag	0.0005-3.5
Al	0.001-11.5
Be	0.0001-0.01
Ca	0.0005-0.3
Cd	0.001-0.07
Cu	0.001-0.2
Fe	0.001-0.05
La	0.002-0.08
Mg	Ref
Mn	0.001-2.4
Ni	0.001-0.05
Pb	0.002-0.1
Si	0.001-0.3
Sn	0.001-0.1
Ti	0.001-0.01
Zn	0.001-7
Zr	0.001-0.7
Li	0.01-11
Y	0.01-5

LIMITS OF DETECTION FOR THE GNR MINILAB 300 in MAGNESIUM BASED MATERIALS