

GNR ANALYTICAL INSTRUMENTS



SUPPLY · SERVICE · SOFTWARE

GNR ANALYTICAL INSTRUMENTS UK - web: www.gnranalytical.co.uk - email: info@gnranalytical.com 0115 937 6686

	NI 01 MINILAB 300 Ni Orientation
Al	0.002-7
B	0.001-3
C	0.002-1
Co	0.005-25
Cr	0.002-35
Cu	0.002-35
Fe	0.002-50
Hf	0.002-1.8
Mg	0.0005-0.1
Mn	0.002-8
Mo	0.002-30
Nb	0.002-6
Ni	Ref
P	0.002-0.1
Pb	0.005-0.1
S	0.002-0.1
Si	0.002-7
Sn	0.002-1
Ti	0.002-6
V	0.002-1
W	0.005-12
Zr	0.002-0.1

LIMITS OF DETECTION FOR THE GNR MINILAB 300 in NICKEL BASED MATERIALS