

GNR ANALYTICAL INSTRUMENTS



SUPPLY · SERVICE · SOFTWARE

GNR ANALYTICAL INSTRUMENTS UK - web: www.gnranalytical.co.uk - email: info@gnranalytical.com 0115 937 6686

TI 01 MINILAB 300 Ti Orientation	
Al	0.001-8
C	0.005-0.1
Cr	0.001-1.2
Cu	0.001-2.5
Fe	0.001-2.5
Mn	0.002-4
Mo	0.001-6
Ni	0.001-0.15
Si	0.001-0.5
Sn	0.002-3
Ti	Ref
V	0.001-6
Zr	0.002-5

LIMITS OF DETECTION FOR THE GNR MINILAB 300 in TITANIUM BASED MATERIALS